

CHEMIA 100

Auto Chemistry Analyzer



Technical Specifications

Operation system

Computer: Built-in operation system - Linux

Screen: 10.4 inch color touch screen

Data storage: 100,000 results

Printer: Built-in and external printer

Others: Soft keyboard

Compatible with keyboard and mouse

Interface:

Ethernet, USB, RS232, VGA, HDMI,

audio, etc.

Working condition

Power supply:

AC 100~240V, 50/60Hz, ≤ 150VA

Temperature: 10-30°C

Humidity: ≤85%

Dimension:

384mm*640mm*410mm (L*W*H)

Weight: 30 Kg

Water consumption: ≤ 1 L/H

Reagent/Sample handling

Reagent/Sample tray:

40 reagent positions

20 sample positions

Reagent volume: 10-300 µL, step by 0.1 µL

Sample volume: 2-70 µL, step by 0.1 µL

Reagent/Sample probe:

Liquid level detection,

vertical & horizontal collision

detection and reagent inventory detection

Probe cleaning:

Auto interior and exterior wash

Optical system

Light source: Halogen tungsten lamp

Photometer: Maintenance-free photometer,

Post spectrophotometry by filters

8 Wavelengths: 340nm, 405nm, 450nm, 510nm,

546nm, 578nm, 630nm, 660nm

Absorbance: 0-4.5 Abs

Hepatic Panel		Immune Panel	
ALB	Albumin	C3	Complement C3
TP	Total Protein	C4	Complement C4
AST	Aspartate aminotransferase	IgA	Immunoglobulin A
ALP	Alkaline phosphatase	IgG	Immunoglobulin G
ALT	Alanine aminotransferase	IgM	Immunoglobulin M
TBA	Total bile acids	RF	Rheumatoid Factor
PA	Prealbumin	ASO	Ascorbic Acid Oxidase
GGT	γ-Glutamyl transferase	Cardiac Panel	
D-BIL	Direct bilirubin	LDH	Lactate dehydrogenase
T-BIL	Total Bilirubin	α-HBDH	α-Hydroxybutyrate dehydrogenase
CHE	Cholinesterase	CK	Creatine Kinase
ADA	Adenosine deaminase	CK-MB	Creatine Kinase Isoenzyme
5'-NT	5'-nucleotidase	MYO	Myoglobin
AFU	α-L-fucosidase	HCY	Homocysteine
Inorganic Panel		Renal Panel	
Ca ²⁺	Calcium	CREA(OX)	Creatinine
Na ⁺	Sodium	UA	Uric Acid
K ⁺	Potassium	UREA	Urea
Cl ⁻	Chlorine	CO ₂	Carbon Dioxide
Mg ²⁺	Magnesium	β ₂ -MG	β ₂ microglobulin
Cu	Cuprum	CysC	Cystatin C
Zn	Zinc	MALB	Microalbuminuria
P	Phosphorus	RBP	Human Retinol Binding Protein
Anemia Panel		Inflammation Panel	
Fe	Iron	CRP	C Reactive Protein
UIBC	Unsaturated iron Binding Capacity	HS-CRP	Hypersensitive-c-reactive-protein
Ferritin	Ferri-albumin	D-Dimer	D-Dimer
TRF	Transferrin	SAA	Serum amyloid A protein
Lipid Panel		Diabetes Panel	
TG	Triglyceride	GLU (OX)	Glucose
CHOL	Total Cholesterol	HBA1C	Glycated Hemoglobin
HDL-C	HDL-Cholesterol	G-6-PD	Glucose-6-Phosphate
LDL-C	LDL-Cholesterol	β-HB	β-Hydroxybutyrate
APOA1	Apolipoprotein A1	FUN	Fructosamine
APOB	Apolipoprotein B	Pancreatitis Panel	
Lpa	Lipoprotein (a)	α-AMY	α-Amylase



Sicuro Bioengineering Ltd


Address: Sicuro House, 19 East West Rd,

Portharcourt, Nigeria.


Email: info@sicurobio.com



CHEMIA 100



Details



Features

Convenient

- Built-in 10.4-inch color touch screen
- Built-in printer
- Solid heating, maintenance-free
- Removable reagent disk

Precision

- Reagent inventory monitoring
- Auto reagent refrigeration
- Accurate probe pipetting technique
- Auto probe washing reduces carry-over
- Simultaneous dual-wavelength measuring...

Economical

- Reaction volume: only 150 μL
- Lower water consumption: $\leq 1\text{L/H}$
- Halogen lamp up to 2000H
- Reusable cuvette
- Fewer consumable

Integrated

- Small footprint on the desk
- Large reagent and sample compartment
- Built-in computer with touch screen
- Complete solution with reagents
- Powerful data transmission with USB, Ethernet, RS232...

Introduction

CHEMIA 100 is a “Mini-size”, fully-auto bench-top chemistry analyzer, the throughput up to 120 tests per hour. It is perfect for low-volume laboratories or as a backup analyzer for mid-to-high volume laboratories.

- » The low cost of ownership of CHEMIA 100 auto chemistry analyzers help laboratories achieve their goals and manage resources effectively without compromising quality.
- » Standardize processes. CHEMIA 100 auto chemistry analyzers features standardized test menus, assay protocols, instrument processes and reference ranges that reduce operator variation and provide consistent quality results.



» System function

General	Benchtop, automatic, random, access, STAT priority
Throughput	Up to 120 tests per hour
Principle	Colorimetry and turbidimetry
Methodology	End-point, fixed-time, kinetic, etc.
Programming	Open or closed system on demand

» Reaction system

Reaction tray	81 reusable cuvettes, auto cuvette washing station
Reaction volume	150 - 750 μL
Reaction temperature	$37 \pm 0.1^\circ\text{C}$
Heating method	Solid heating

» Calibration and control

Calibration	K-factor, Linear, Spline, Logit-Log 4P, Logit-Log 5P, Exponential, Polynomial
Control	Westgard multi-rule, Cumulative sum, check, Twin plot, L-J graph