

# IMMUNIA-100

Immunofluorescence Analyzer



- Rapid Testing
- Easy to operate
- Reliable accuracy
- Complete test item

**SICURO**

# Immunia-100

## Immunofluorescence Analyzer

■ Immunia-100 Immunofluorescence Analyzer is an advanced POCT analyzer, intended to deliver accurate testing results rapidly.



### ■ Application

Immunia-100 is an essential analyzer for Emergency Department, ICU, NICU, Cardiovascular Department or anywhere treatment and research needs at the point of care.



### ■ Specification

**Weight**  
4kg

**Work temperature**  
10~30°C

**Interface**  
5" touch screen  
Built-in thermal printer  
USB port / Ethernet  
LIS/HIS compatible

**Memory**  
50,000 results


**Power supply**  
100-240V AC, 50~60 Hz

**Dimension**  
280mm\*260mm\*140mm (L\*W\*H)

## Advantages

### Rapid Testing

Delivering results within 3-15 mins



### Easy to operate

Touch Screen operation, Barcode scanner is optional



### Reliable accuracy

Multi-points correlation, the testing results are reliable



### Complete test item

Inflammation and cardiac markers were tested, including H-FABP



- Rapid screening, results available in just 5 minutes
  - Full color screen
  - Regular QC testing
  - Supports multiple languages
  - The reserve capacity exceeds 10,000 servings
  - Printing capability
-

## Testing Menu

### Thyroid

TSH T3 T4 FT3 FT4

### Cardiac

NT-proBNP D-Dimer cTnl Myoglobin BNP CK-MB cTnI/CK-MB/Myo CTnl / NT-proBNP / D-Dimer

### Inflammation

PCT IL-6 Hs-CRP SAA PCT & IL-6 Hs-CRP CRP & Hs-CRP ( 2 in 1) SAA & W-CRP ( 2 in 1)

### Fertility

AMH  $\beta$ -HCG LH FSH Prolactin Progesterone Testosterone

### Renal Panel

CYS-C NGAL MALB Cortisol

### Tumor Markers

AFP CEA f-PSA T-PSA CA125 CA15-3 CA19-9

### Anemia

Ferritin Folate

### Specific protein

IgE

### Nerve and brain injury

S100- $\beta$

### Diabetes

HbA1c

### Bone Metabolism

Vitamin D

### Drug Abuse

MET<sup>+</sup> KET<sup>+</sup> MOP<sup>+</sup> MDMA<sup>+</sup> COC<sup>+</sup> MET<sup>+</sup>/KET<sup>+</sup>/MOP<sup>+</sup> THC<sup>+</sup> K2<sup>+</sup>

## Sicuro Bioengineering Ltd

### Web

<http://www.sicurobio.com>

### Email

[info@sicurobio.com](mailto:info@sicurobio.com)

### Add

Sicuro House, 19 East West Rd, Portharcourt, Nigeria.