

HBSAG RAPID TEST

BASED ON A HIGHLY SENSITIVE TWO-SITE SANDWICH IMMUNOASSAY

The complex antigen found on the surface of Hepatitis B virus (HBV) is called HBsAg. The presence of HBsAg in whole blood, serum or plasma is an indication of an active Hepatitis B infection, either acute or chronic. In a typical Hepatitis B infection, HBsAg will be detected 2 to 4 weeks before the ALT level becomes abnormal and 3 to 5 weeks before symptoms or jaundice develop. HBsAg has four principal subtypes: adw, ayw, adr and ayr. Because of antigenic heterogeneity of the determinant, there are 10 major serotypes of Hepatitis B virus.

Ordering Information

Cat. No.: SIC-RDT-HBSG-202022

Product: HBsAg Rapid Test

Specimen: WB/S/P

Format: Dipstick

Kit Size: 50 T

Product Advantages

- ✓ Fast results in **15-30** minutes
- ✓ Easy visual interpretation
- ✓ High accuracy

Convenient Procedure

