

MULTI 12 DRUGS RAPID TEST 1-STEP CUP

COMPREHENSIVE 12-PANEL SCREENING IN ONE STEP

Analytes	Cut off	Analytes	Cut off
Methadone (MTD)	300 ng/mL	Marijuana (THC)	50 ng/mL
Morphine (MOP)	300 ng/mL	Amphetamine (AMP)	1000 ng/mL
Methamphetamine (MET)	1000 ng/mL	Tramadol (TML)	100 ng/mL
Ketamine (KET)	1000 ng/mL	Diazepam (DIA)	300 ng/mL
Ecstasy (MDMA)	1000 ng/mL	Cannabinol (CNB)	500 ng/mL
Cocaine (COC)	300 ng/mL	Clonazepam (CLO)	400 ng/mL

Ordering Information

Cat.no.: SIC-DOA12-01
 Product Description:
 Multi 12 Drugs Rapid Test 1-Step Cup
 Specimen: Urine
 Format: Cup
 Pack: 10 T
 CE Status: CE

Product Advantages

- ✓ Up to 12 analytes in one go
- ✓ High Accuracy
- ✓ Fast result in 5 minutes
- ✓ Convenient 1- step procedure

Convenient Procedure



Read the drug strips in 5 minutes

