ELISURE

CLOSTRIDIUM DIFFICILE GDH RAPID TEST

TARGET THE GDH ENZYME: HIGH-SENSITIVITY ANTIGEN CAPTURE

Toxinogenic strains of *Clostridium difficile* cause infections from mild diarrhea to pseudomembranous colitis, potentially leading to death. Disease is caused by two toxins produced by toxinogenic strains of C.difficile: Toxin A (tissue-damaging enterotoxin) and Toxin B (cytotoxin). The use of Glutamate Dehydrogenase (GDH) as an antigen marker of C.difficile proliferation has been shown to be very effective because all strains produce high amount of this enzyme.
Clostridium difficile GDH Rapid Test allows the specific detection of *C. difficile's* GDH in stool specimen. Samples with a positive result should be investigated further to test for toxigenicity of the bacteria.


Ordering Information

Cat. No.: SIC-RDT-CDIF-202064
Product: *Clostridium difficile* GDH Rapid Test
Specimen: Feces
Format: Cassette
Kit Size: 10 T

Product Advantages

- ✓ Fast results in **10** minutes
- ✓ Easy visually interpretation
- ✓ High Accuracy

Convenient Procedure





Sicuro Bioengineering Ltd.
Address: 19 East West Road, Port Harcourt - Nigeria
Email: info@sicurobio.com
Website: www.sicurobio.com

ISO ISO

* CE Marking in progress