ELISURE

FOB RAPID TEST

SCREEN FOR COLORECTAL CANCER

Many diseases can cause hidden blood in the feces. This is also known as Fecal Occult Blood (FOB), Human Occult Blood, or Human Hemoglobin. In the early stages, gastrointestinal problems such as colon cancer, ulcers, polyps, colitis, diverticulitis, and fissures may not show any visible symptoms, only occult blood.

The FOB Rapid Test is a rapid test to qualitatively detect low levels of Fecal Occult Blood. The test uses a double antibody sandwich assay to selectively detect Fecal Occult Blood at 50 ng/mL or higher, or 6 µg/g feces. In addition, unlike guaiac assays, the accuracy of the test is not affected by the diet of the patients.


Ordering Information

Cat. No.: SIC-RDT-FOB-202037
Product: FOB Rapid Test
Specimen: Feces
Format: Cassette
Kit Size: 25 T

Product Advantages

- ✓ Fast results in **5** minutes
- ✓ Easy visually interpretation
- ✓ High accuracy

Convenient Procedure





Sicuro Bioengineering Ltd.
Address: 19 East West Road, Port Harcourt - Nigeria
Email: info@sicurobio.com
Website: www.sicurobio.com

ISO ISO

* CE Marking in progress