ELISURE

TETANUS AB RAPID TEST

CHECK THE IMMUNITY AGAINST LIFE--THREATENING TETANUS IN TIME

Tetanus is caused by an infection with the bacterium Clostridium tetani, which is commonly found in soil, saliva, dust, and manure.
The bacteria generally enter through a break in the skin such as a cut or puncture wound by a contaminated object.
Tetanus occurs in all parts of the world , as C. tetani spores are ubiquitous, but is most frequent in hot and wet climates.
Typical Symptoms of Tetanus : Muscle spasms, Fever, Headache
Testing with Tetanus rapid tests lets the doctors understand the immunity status and decide on vaccination accordingly.
The Tetanus Rapid Test Cassette (Whole Blood/Serum/Plasma) is a rapid chromatographic immunoassay for the qualitative detection of antibodies to Tetanus toxin in whole blood, serum, or plasma to aid in the diagnosis of Tetanus toxin infection.

Ordering Information

Cat. No.: SIC-RDT-TET-202070

Product: Tetanus Rapid Test

Specimen: WB/S/P


Format: Cassette

Kit Size: 25 T

Product Advantages

- ✓ Fast results in **10** minutes
- ✓ Simple Procedure
- ✓ High accuracy
  (Sensitivity: 94.1% ; Specificity: 98.9% ; Accuracy: 97.8%)

Convenient Procedure





Sicuro Bioengineering Ltd.
Address: 19 East West Road, Port Harcourt - Nigeria
Email: info@sicurobio.com
Website: www.sicurobio.com

ISO ISO

* CE Marking in progress