ELISURE

TYPHOID IGG/IGM RAPID TEST

DISTINGUISHING ACUTE INFECTION FROM PAST EXPOSURE

Typhoid fever is caused by S. typhi, a Gram-negative bacterium. World-wide an estimated 17 million cases and 600,000 associated deaths occur annually. Patients who are infected with HIV are at significantly increased risk of clinical infection with S. typhi. 1-5% of patients become chronic carrier harboring S. typhi in the gallbladder.

The Typhoid IgG/IgM Rapid Test is a simple and rapid laboratory test. The test simultaneously detects and differentiates the IgG and the IgM antibodies to S. typhi specific antigen in whole blood, serum or plasma thus aid in the determination of current or previous exposure the S. typhi.


Ordering Information

Cat. No.: SIC-RDT-TYP-202094
Product: Typhoid IgG/IgM Rapid Test
Specimen: WB/S/P
Format: Cassette
Kit Size: 25 T

Product Advantages

- Fast results in **15** minutes
- Easy visually interpretation
- High accuracy

Convenient Procedure



Sicuro Bioengineering Ltd.
Address: 19 East West Road, Port Harcourt - Nigeria
Email: info@sicurobio.com
Website: www.sicurobio.com

ISO ISO

* CE Marking in progress