ELISURE

VIBRIO CHOLERAE O1 (VC O1) RAPID TEST

AID IN THE DIAGNOSIS OF VIBRIO CHOLERAE INFECTION

Cholera is an important public health problem in many parts of Asia, Africa and Latin America. Globally, 3–5 million cases and over 100,000 deaths occur annually due to cholera. Countries facing complex emergencies are more vulnerable to cholera outbreaks. The etiological agent of cholera has been identified as Vibrio cholerae (V. cholerae), a gram negative bacterium, which is generally transmitted to humans via contaminated water and food. The Vibrio cholerae O1 (VC O1) Rapid Test is a rapid chromatographic immunoassay for the qualitative detection of Vibrio cholerae O1 in human feces, providing results in 10 minutes.


Ordering Information

Cat. No.: SIC-RDT-VC-202065
Product: *Vibrio cholerae O1 (VC O1)* Rapid Test
Specimen: Feces
Format: Cassette
Kit Size: 10 T

Product Advantages

- ✓ Fast results in **10** minutes
- ✓ Easy visually interpretation
- ✓ High accuracy

Convenient Procedure





Sicuro Bioengineering Ltd.
Address: 19 East West Road, Port Harcourt - Nigeria
Email: info@sicurobio.com
Website: www.sicurobio.com

ISO ISO

* CE Marking in progress